DE

EN

FR

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
THANK YOU FOR BUYING A PAIR OF HIGH-QUALITY FST SPRING SCISSORS.

Every single one of our instruments has been tested by our surgical technicians and adapted to your needs at the micrometre level.

The result: extremely fine cutting edges, joints and screws – all of which require the same level of care when it comes to handling and maintenance as they did when they were being produced.

- You should clean, rinse and dry the instruments carefully after use.
- For cleaning, we recommend a soft toothbrush, always brushing towards the tip.
- After every other sterilisation, you should wet the screws, joints and moving parts of the scissors with a drop of our instrument oil FST # 29055-00 or with our oil pen FST # 29055-12. Please do not use WD-40 or other common industrial oils for this.

This way, you can avoid corrosion and damage to the fine springs, rongeurs, scissors, needle holders and clamps – and ensures that you can always 100 % rely on your instruments.

If you have any other technical questions, please don't hesitate to contact us.

STERILISATION TIPS

DE

EN

FR

STEAM STERILISATION / AUTOCLAVE PROCESS

The autoclave process uses pressurised moist heat for sterilisation. It's important that the steam used comes from sterilised water and is not contaminated.

STERILISATION TIPS – CONTINUED

- You can sterilise your instruments individually or as a set.
- Individual instruments:
Disposable cases or sterilisation bags are ideal for this (FST 28000-07, 28000-10, 28000-12). Make sure that the bag is sufficiently large – the instrument should be in an open (unlocked) position.
- Instrument sets:
The instruments should all be open / unlocked here too. If it is necessary to have two layers, the heavy instruments should be at the bottom.

IMPORTANT: Please DO NOT mix instruments of different or unknown alloys in the autoclave.

- Whether you're sterilising your instruments individually or as a set, do not close / lock your instruments during autoclaving. Otherwise, the steam won't be able to reach and sterilise the adjoining metallic surfaces. In addition, thermal expansion during autoclaving may result in cracks in the joint areas.

STERILISATION TIPS – CONTINUED

- Do not overload the autoclave chamber. Doing so would impede steam penetration, meaning the instruments may not be 100 % sterilised.
- Tip: Place a towel beneath the pan to catch excess moisture during autoclaving. Make sure that the towel has no washing detergent residue on it, and is pH-neutral when submerged in water.
- At the end of the autoclaving cycle (before the drying cycle), unlock the autoclave door and open it slightly. Then allow the drying cycle to run for the period of time recommended by the autoclave manufacturer.

PLEASE NOTE: If you open the autoclave door fully before the drying cycle, cold air will flow into the chamber – and this may lead to condensation forming on the instruments and cause water spots on the instruments.

| STERILISATION TIPS – CONTINUED

ETHYLENE OXIDE STERILISATION (ETO STERILISATION)

Some instruments, such as catheters, cannulas, etc., cannot withstand high temperatures or high humidity levels. Ethylene oxide sterilisation (ETO sterilisation) is recommended for these instruments. This is a process composed of four primary variables: gas concentration, humidity, temperature and time. You will find suitable sterilization accessories in our catalog under the heading Instrument Care.

OXIDATION (RUST)

FST provides high quality, stainless steel surgical instruments that are resistant to oxidation or rust. Each individual instrument is manufactured in strict accordance with European standards for stainless steel medical goods. You should maintain these instruments with great care in order to extend their service life. Stainless steel is, of course, corrosion-resistant. Nevertheless, improper handling may result in rusting or discolouration – if stainless steel is subjected to physical, thermal or chemical damage, it may lose some of its rust-proof properties.

REPAIR TIPS


Our forceps are highly-precise instruments that have been hand-crafted to meet your needs by experienced surgical technicians – so repair work may require a little practice.

We recommend the following procedure: Try out the following steps for the first time on old forceps that you don't need anymore.


STRAIGHTENING FORCEPS:

- Using pliers, grip the forceps right at the spot where they are bent.
- Apply more pressure to the tip on the inside (the surface that is touching the other tip) and gently squeeze the pliers.
- Please note that stainless steel and carbon steel are harder alloys and are more likely to break at the bent points than become straight again. Dumoxel, Dumostar and titanium are more flexible alloys.


ALIGNING FORCEPS:

- First, determine where the points of the forceps diverge.
- Then use the same technique on the pliers as in the straightening instructions above, but with more pressure on the side where they are bent – exactly opposite the tip.
- Bend the tip further inwards. The forceps will now close properly again.
- If the tips are not the same length, keep the forceps closed and run the tips gently over a sharpening stone until the protruding end has been removed.


SHARPENING FORCEPS:

- Use the grooved sharpening stone (FST 29008-22) for sharpening forceps.
- Use the supplied O ring to keep the tips of the forceps together.
- Hold the forceps gently by the edges and move them gently back and forth on the sharpening stone.
- Always use the mineral oil supplied (FST 29055-00) when using the sharpening stone.
- Continuously observe the progress of your work under a microscope – and make sure that you sharpen both sides of the edge.

SHARPENING FORCEPS:

- Keep the forceps as parallel to the sharpening stone as possible and repeat the sharpening process on both sides – the more parallel you work, the less abrasion there will be.
- Continuously change the side of the sharpening stone you are using, and clean the dust from the stone with mineral oil every so often.
- Use the grooves in the sharpening stone for the final fine sharpening. To maintain equal tips, rotate the instrument back and forth as you move it over the stone.

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Fine Science Tools (USA) Inc.
Fine Science Tools (CAN) Inc.
Fine Science Tools GmbH